


Rev 2.0

# **ResMed Ventilators and COVID-19**

## *Information on applications in the treatment of patients with COVID-19*

The information contained in this document is current as of March 23, 2020, and is based on currently available information that will continue to change over time. The information in this guide with respect to treatment is believed to have a reasonable basis. ResMed assumes no obligation to update the information in this presentation, whether as a result of new information or future events.



## Table of Contents

Executive Summary.....	3
Use of Ventilation for Patients with COVID-19 .....	3
Invasive ventilation is important for critically ill patients .....	3
Non-invasive ventilators can help in decompressing hospital census.....	4
Utilizing oxygen.....	5
Protecting healthcare workers to mitigate infection risk.....	5
Clinical training for non-invasive ventilation.....	6
ResMed Ventilators and Bi-level Devices – Application and Features .....	7
What is the difference between non-invasive and invasive ventilation?.....	8
What is a bi-level device? .....	8
Invasive and non-invasive ventilators.....	9
Astral 100/150 .....	9
Stellar 100/150 .....	9
Non-invasive ventilators.....	9
GA ST series .....	9
Bi-level devices for non-invasive ventilation .....	10
Lumis and AirCurve.....	10
Flexo ST series .....	10
Summary of device specifications .....	10
Reprocessing of ResMed Devices.....	11
References and Notes.....	12
APPENDIX – Device Specifications.....	13



## Executive Summary

As a medical technology manufacturer of respiratory devices, including ventilators, ResMed has mobilized a global task force in response to the COVID-19 pandemic. Our guiding principle in responding to the COVID-19 pandemic is the preservation of life throughout the globe.

Access to ventilators will be important during this pandemic and ResMed has taken the necessary steps to prioritize manufacturing resources for devices that support higher acuity patients. However, as ventilator resources become stretched due to high demand, alternative ventilation options such as bi-level devices will be crucial to stabilizing or sustaining patients who require respiratory support.

This document aims to assist healthcare officials in understanding the application of ResMed's devices in providing ventilation support to patients with clinical syndromes due to COVID-19 infection. It is based on current information, which is changing rapidly. **The purpose is not to direct clinical practice but to provide clear information on the available ResMed products and their application.**

## Use of Ventilation for Patients with COVID-19

Patients affected by COVID-19 have a range of symptoms and varying degrees of severity of illness. All patients should be assessed thoroughly, with the decision to monitor a particular patient in the inpatient or outpatient setting depending on clinical presentation and severity of illness, the patient's ability to engage in monitoring, home isolation, and the risk of transmission.<sup>1</sup>

According to the Center for Disease Control (CDC), patients with mild illness on presentation may not initially require hospitalization; however, their symptoms may worsen with progression to lower respiratory tract disease.<sup>1</sup> Many are diagnosed with pneumonia and subsequently develop hypoxemic respiratory failure or acute respiratory distress syndrome (ARDS). These patients are primarily managed with mechanical ventilation, both invasive and non-invasive, often with supplemental oxygen.<sup>2</sup> Published information shows that 6% of patients with COVID-19 require ventilation,<sup>3</sup> with numbers rising to 89% for those in the intensive care unit (ICU).<sup>4</sup> Of these, 47.2% received invasive ventilation and 41.7% received non-invasive ventilation.<sup>4</sup>

## Invasive ventilation is important for critically ill patients

Several guidance documents for the management of critically ill patients with COVID-19 have been published;<sup>1,2</sup> these are by and large based on the usual management of viral pneumonia with respiratory failure, with additional precautions to reduce risk of transmission. Current guidance from the World Health Organization (WHO) and others recommend early consideration for invasive ventilation for patients with severe COVID-19 who develop ARDS,<sup>2,5</sup> and a strong preference for early use of invasive ventilation over non-invasive ventilation where appropriate and possible.<sup>6</sup> Additional information from the Handbook of COVID-19 Prevention and Treatment indicated that some severe patients progress to ARDS rapidly and intubation should be performed as early as possible if improvement in respiratory distress symptoms or PaO<sub>2</sub>/FiO<sub>2</sub> is not observed.<sup>7</sup> Evidently, access to invasive ventilators will be critical



during this pandemic. However, the dramatic increase in patients requiring ventilator therapy may lead to ventilator scarcity. Emerging data from China show that only 25% of patients who died had received invasive ventilation, suggesting that ventilation resources may not have been available for critically ill patients.<sup>8</sup>

Governments around the world, recognizing the impending shortage, are working hard to prepare and deploy resources as quickly as possible. The UK government released a Specification on 22 March 2020 for minimally clinically acceptable ventilator systems that are considered suitable for hospital use in treating COVID-19. This specification guidance outlines the need for volume controlled (Pressure Regulated Volume control or pressure controlled ventilation) ventilation as well as bi-level devices while highlighting the monitoring and safety features expected of ventilators used in critical care situations.<sup>9</sup>

The U.S. Food and Drug Administration (FDA) released an enforcement guidance on 22 March 2020 which provides a policy to help expand the availability of ventilators and other respiratory devices during this pandemic. The guidance describes the FDA's intention to exercise enforcement discretion for certain deviations, such as the use of ventilators outside their cleared environment of use, and the use of devices indicated for sleep apnea (including noncontinuous ventilators delivering CPAP or bi-level positive airway pressure) to treat patients with respiratory insufficiency, provided that appropriate design mitigations are in place to minimize aerosolization.<sup>10</sup>

## **Non-invasive ventilators can help in decompressing hospital census**

In a recent discussion paper, the National Academy of Medicine stated that the use of bi-level devices could be a way to forestall the need for intubation and reduce days on a ventilator.<sup>11</sup> Bi-level devices provide non-invasive ventilation only, and could be used to stabilize less severe patients while allowing allocation of ventilators with invasive capabilities to critical cases, or sustain a patient long enough for an invasive ventilator to become available.<sup>12</sup>

A Philips Healthcare webinar on respiratory management for severe COVID-19 patients showed the impact non-invasive ventilation devices have had in treating patients with COVID-19 in China and Italy.<sup>13</sup> During the webinar, Professor Cheng Rongchang discussed cases of COVID-19 patients using non-invasive ventilation with bi-level devices; he shared that *“four out of sixteen patients in the ICU were managed with non-invasive ventilation and avoided intubation”*. Professor Stefano Nava, the chief of the Respiratory and Critical Care Unit at Sant' Orsola Hospital in Bologna, Italy, shared his experiences of desperate shortages of ventilators and the need to triage patients by their respiratory symptoms, with non-invasive ventilation used for those who needed more than just supplementary oxygen and moving to intubation if non-invasive ventilation could not improve their oxygen levels. These experiences of physicians treating severe and critical COVID-19 patients in China and Italy has confirmed the ability of triaging some patients to non-invasive ventilation. As always, patient selection is important and all patients must be monitored closely for deterioration.

While non-invasive ventilation is not ideal for the most severe or critical cases of COVID-19, this therapy type could be useful in decompressing hospital census. Based on current literature and anecdotal



evidence from the field, it appears that bi-level devices or non-invasive ventilators may be useful in the following scenarios:

1. When a patient needs support for respiratory insufficiency, but has not deteriorated into more severe hypoxemia, ARDS, or any other clinical scenario where invasive ventilation is more appropriate.
2. To decompress ICU census, after a patient is extubated and recovering from invasive ventilation, and help mitigate the risk of reintubation, to allow an invasive ventilator to be cleaned and serviced and circulated back into use on another patient.
3. To decompress hospital census, allowing patients who still need some respiratory support to transition to the home or other non-hospital facilities (e.g. rehabilitation center, skilled nursing facility, etc.) on devices that are already used in millions of homes every day.
4. The patient needs intubation but there are no more invasive ventilators, non-invasive ventilation could be useful to bridge the time to intubation. This situation is the highest risk and not ideal therefore the patient must be closely monitored. Some evidence suggests this may only be helpful for a short time (e.g. less than 2 hours).<sup>7</sup>

In the Philips Healthcare webinar, Professor Stefano Nava recounted the experience from Italy where ICU beds and equipment were fully utilized and entire hospitals were dedicated for the treatment of patients with COVID-19.<sup>13</sup> Practicality became the key to dealing with the influx of patients requiring treatment, with those with mild symptoms being sent home for monitoring.<sup>13</sup> For countries facing this public health emergency, and where ICU capacity is insufficient to meet demand, it will be important to provide ventilation solutions to support patients into the sub-acute or out of hospital care environments.

## Utilizing oxygen

Hypoxemia can present due to impaired respiratory functions by COVID-19. Oxygen supplementation treatment may correct hypoxemia, relieving secondary organ damage caused by respiratory distress and hypoxemia. Evidence suggests that controlled oxygen therapy may be helpful,<sup>5</sup> and this can be done with invasive or non-invasive ventilation by oxygen entrainment into the circuit or patient interface.


## Protecting healthcare workers to mitigate infection risk

Some concerns have emerged regarding the risk of dispersion of aerosolized virus when utilizing non-invasive ventilation. However, evidence suggests that non-invasive ventilation procedures are more likely to produce large droplets (>10  $\mu\text{m}$ ) rather than aerosols, and that these are largely confined to within one meter due to their large mass.<sup>14</sup> This suggests that the risk of droplet dispersion as a result of use of non-invasive ventilation or bi-level devices may not be that different to that of any COVID-19 patient in the hospital who is coughing or sneezing.

Additionally, an experts' panel determined that non-invasive ventilation systems with a good interface fitting do not create widespread dispersion of exhaled air.<sup>15</sup> Recommendations have been published to support good mask fit to reduce aerosols, including use of full face masks.<sup>16</sup> Nonetheless, the risk of

aerosol dispersion needs to be mitigated with appropriate isolation of patients and the use of Personal Protective Equipment (PPE) for healthcare workers, such as N95 masks/respirators and eye protection,<sup>5</sup> which are standard protective equipment in a COVID-19 ICU.<sup>17</sup>

Recommended methods for reducing risk of virus spread include using suitable masks and filters, appropriate PPE and isolation techniques (see **Figure 1**).



**Figure 1. Mitigating risk of droplet dispersion when using non-invasive ventilation.**<sup>1, 7</sup>

## Clinical training for non-invasive ventilation

The principles of invasive ventilation and non-invasive ventilation are very similar, so training healthcare professionals (e.g. anesthesiologists, emergency physicians, intensivists, nurses, and respiratory therapists) who are well-versed in invasive ventilation to provide non-invasive ventilation should create few burdens. The two primary changes between invasive ventilation and non-invasive ventilation are the circuit configurations and the modes of ventilation utilized. There are fewer required circuit components when initiating non-invasive ventilation, so the burden is expected to be less than that of initiating invasive ventilation.


The most challenging part of setting up non-invasive ventilation is training care givers to properly fit a mask. Poor mask fit can cause discomfort or intolerance, spread exposure to health care providers and reduce therapy effectiveness. The ability to anticipate, prevent and manage mask-related problems will be important for non-invasive ventilation success.

ResMed offers online tutorials of device setups, quick setup guides, suggested non-invasive ventilation settings and remote webinar trainings for users, all tools which largely allow for independent clinician setup of devices.

The method of ventilation is an important clinical decision to be made by the treatment team under rapidly evolving clinical guidelines for COVID-19 patients, availability of ventilation technology, clinical settling and availability of personal protective equipment for healthcare workers. The following information about ResMed devices is designed to assist and inform these decisions by clarifying applications and features of different devices.

## ResMed Ventilators and Bi-level Devices – Application and Features

ResMed manufactures a range of ventilators and bi-level devices. These devices are indicated for hospital and home use, and have the flexibility for use in various clinical scenarios (see **Figure 2**). It should be noted, however, that these are **not** the same as ventilator equipment typically used in high acuity situations in hospital ICUs.



**Figure 2. ResMed and Curative ventilators and bi-level devices in various clinical scenarios**

The FDA in their COVID-19 ventilator enforcement guidance encourages applications to distribute in the U.S. ventilators and respiratory devices that have been approved in other jurisdictions;<sup>10</sup> this increases the ability of manufacturers like ResMed to respond.

The ResMed range of ventilators and bi-level devices capable of providing respiratory support to patients at various stages of dependency is shown below (see **Figure 3**).




Figure 3. ResMed and Curative ventilators and bi-level devices

## What is the difference between non-invasive and invasive ventilation?

**Non-invasive ventilation (NIV)** is a form of mechanical ventilation where air is delivered to the patient through a mask or mouthpiece. **Invasive ventilation** is used when sufficient ventilation cannot be achieved using non-invasive methods so air is delivered through a tube inserted into the trachea either via intubation or tracheotomy.

## What is a bi-level device?

Bi-level devices are considered ventilatory equipment and provide non-invasive ventilation only. These devices deliver two distinct pressures, one for inhalation and the second for exhalation; the change in pressures leads to flow of air in and out of the lungs. Bi-level devices are not cleared or approved by local regulatory authorities for use with intubation or tracheotomy. In most European and North American countries, most of these devices are not utilized by physicians for invasive ventilation because they lack alarms, sophisticated monitoring or modes, and their pressure generation is not sufficient to care for those with significant lung injury. However, some advanced bilevels are suitable for invasive ventilation and cleared in several jurisdictions for that use.



## Invasive and non-invasive ventilators

### Astral 100/150

Of the ResMed range of invasive and non-invasive ventilators, the Astral life-support ventilator provides the most comprehensive set of modes and settings and delivers both pressure and volume ventilation. The ResMed Astral 150 is the most comprehensive—it comes with all the standard features of the standard Astral 100 plus double limb circuit capabilities, which complements the use of an inspiratory and expiratory antibacterial/antiviral filter and a non-vented mask to reduce risk of contamination to healthcare professionals. The device allows addition of supplemental oxygen (no oxygen blender; low pressure oxygen only) at the air inlet up to 30L/min and provides monitoring of FiO<sub>2</sub>. The ability to use an active circuit may help reduce particle spread. The device can generate higher pressures needed to care for an acutely ill patient. The option to have four pre-set settings may make it easier for less experienced staff taxed during the global health crisis to run the device. These devices also include remote monitoring that can facilitate “telehealth” or remote management of patients, as well as the ability to centralize telemonitoring of devices in a “war room” configuration

### Stellar 100/150

The ResMed Stellar 100/150 is a non-invasive ventilator with invasive capabilities when combined with the ResMed leak valve, and is indicated for ventilation of non-dependent, spontaneously breathing patients. The Stellar 100/150 can be safely used in those recovering from acute lung injury such as ARDS or those with milder underlying lung disease. The device also allows addition of supplemental oxygen at the air inlet up to 30L/min and provides monitoring of FiO<sub>2</sub> with an additional FiO<sub>2</sub> monitoring sensor attached. This device has an internal battery and can be used for transport within a hospital. Physicians consider this device a good option for step down units.

## Non-invasive ventilators

### GA ST series

The GA ST series is a non-invasive ventilator indicated for ventilation of spontaneously breathing individuals. It is equipped with two oxygen input methods, high-pressure oxygen and low-pressure oxygen. The device allows addition of supplemental oxygen at the air inlet and is capable of providing up to 100% FiO<sub>2</sub>. The maximum oxygen flow rate with low-pressure is 30 L/min. An oxygen mixer is also available for greater oxygen concentration accuracy. The use and replacement of the air inlet filter and main flow bacteria filter between patients and at regular intervals will reduce the risk of patient or ventilator contamination. This device is sold under the company name Curative, a ResMed family company.



## Bi-level devices for non-invasive ventilation

### Lumis and AirCurve

The Lumis and AirCurve range of ResMed devices<sup>18</sup> are bi-level devices indicated to provide non-invasive ventilation for patients with respiratory insufficiency. A "backup" rate can be set to ensure that patients still receive a minimum number of breaths per minute if they fail to breathe spontaneously. The ST-A variant comes with fixed and adjustable alarms to alert the user/caregiver in case of therapy issues and may be more appropriate for the ward environment than the ST variant, which is possibly more suitable for when a patient is discharged to a home environment. If necessary, supplemental oxygen up to 15 L/min can be connected to the air outlet of the Lumis and AirCurve range, but monitoring of FiO<sub>2</sub> is not done by the device. An additional oximetry adapter can be attached to measure SpO<sub>2</sub> if needed. These devices also include remote monitoring that can facilitate "telehealth" or remote management of patients, as well as the ability to centralize telemonitoring of devices in a "war room" configuration.

### Flexo ST series

The Flexo ST series is a bi-level device that provides non-invasive ventilator support for patients with respiratory insufficiency. Like the Lumis and AirCurve range, a timed backup mode allows the minimum number of breaths per minute to be set in the event that the patient slows or ceases respiratory efforts. The highest end of the Flexo ST series provides inspiratory pressures up to 30 cmH<sub>2</sub>O and includes adjustable alarms to alert the user/caregiver in case of therapy issues. The Flexo device is not designed or intended to be used with supplemental oxygen. This device is sold under the company name Curative, a ResMed family company.

## Summary of device specifications

Key specifications of ResMed ventilators and bi-level devices are provided in the table below (see **Table 1**). Comprehensive information on these devices, their specific indications for use, modes and features is set forth in the Appendices.



**Table 1. Key specifications of ResMed and Curative ventilators and bi-level devices**

	<b>Astral 100/150</b>	<b>Stellar 100/150</b>	<b>GA ST*</b>	<b>Lumis/AirCurve ST-A</b>	<b>Lumis/AirCurve ST</b>	<b>Flexo ST*</b>
<b>Invasive Capability</b>	Yes	Yes (leak valve with tracheostomy)	No	No	No	No
<b>Non-invasive Capability</b>	Yes	Yes	Yes	Yes	Yes	Yes
<b>Max Pressure (cmH<sub>2</sub>O)</b>	50	40	40	30	25	30
<b>O<sub>2</sub> Entry Location</b>	Inlet	Inlet	Inlet	Outlet	Outlet	N.A.
<b>Max O<sub>2</sub> Flow Rate (L/min)</b>	30	30	30	15	15	N.A.
<b>FiO<sub>2</sub> Monitoring/Alerts</b>	Yes	Yes	Yes	No	No	No
<b>Humidification</b>	External	Integrated or External	External	Integrated or External	Integrated or External	External
<b>Alarms</b>	Yes	Yes	Yes	Yes	No	Yes
<b>Internal battery</b>	Yes	Yes	Yes	No	No	No
<b>External battery</b>	Yes	Yes		Yes	Yes	No
<b>Telemonitoring</b>	Yes	Yes	No	Yes	Yes	No
<b>Modes</b>	CPAP, (S)T, P(A)C, (A)CV, P(A)CV, P-SIMV, V-SIMV, PS, iVAPS	CPAP, S, ST (optional iBR), T, PAC, iVAPS	CPAP, S, T, ST, APCV, CPAV, TVV-t, TVV-ST, TVV-APCV	CPAP, S, ST (optional iBR), T, PAC, iVAPS	<i>AirCurve</i> : CPAP, S, ST, T <i>Lumis</i> : CPAP, S, ST (optional iBR), T, PAC, iVAPS	CPAP, S, T, ST, APCV, TVV

\*GA ST and Flexo ST are devices offered by Curative, a subsidiary of ResMed.


## Reprocessing of ResMed Devices

The novel coronavirus that causes the disease COVID-19, SARS-CoV-2, is an enveloped virus. This means that it can be inactivated with the appropriate disinfectant product when used according to the label directions. The U.S. Environmental Protection Agency has published a list of disinfectants that meets its criteria for use against SARS-CoV-2.<sup>19</sup> This authorized list of disinfectants is comprised of many commonly used disinfectants and is being actively updated as more information emerges.

Published guidance from health authorities reinforce the need to maintain standard cleaning and disinfection procedures.<sup>20</sup> For each ResMed device, these cleaning and disinfection procedures are provided in the device’s associated clinical guide, user guide or service manual. To prevent cross-contamination, antibacterial filters are used on air intake and circuits, and circuit accessories are to be replaced or sterilized. Instructions are also provided in the materials and method for cleaning surfaces.

## References and Notes

- <sup>1</sup> Centers for Disease Control and Prevention. Interim Clinical Guidance for Management of Patients with Confirmed Coronavirus Disease (COVID-19). Retrieved from <https://www.cdc.gov/coronavirus/2019-ncov/hcp/clinical-guidance-management-patients.html> March 22, 2020
- <sup>2</sup> WHO. Clinical management of severe acute respiratory infection (SARI) when COVID-19 disease is suspected: interim guidance, 13 March 2020. Retrieved from [https://apps.who.int/iris/bitstream/handle/10665/331446/WHO-2019-nCoV-clinical-2020\\_4-eng.pdf?sequence=1&isAllowed=y](https://apps.who.int/iris/bitstream/handle/10665/331446/WHO-2019-nCoV-clinical-2020_4-eng.pdf?sequence=1&isAllowed=y) March 17, 2020.
- <sup>3</sup> Guan WJ, Ni ZY, Hu Y, et al. Clinical Characteristics of Coronavirus Disease 2019 in China [published online ahead of print, 2020 Feb 28]. *N Engl J Med*. 2020
- <sup>4</sup> Wang D, Hu B, Hu C, et al. Clinical Characteristics of 138 Hospitalized Patients With 2019 Novel Coronavirus-Infected Pneumonia in Wuhan, China [published online ahead of print, 2020 Feb 7]. *JAMA*. 2020
- <sup>5</sup> Murthy S, Gomersall CD, Fowler RA. Care for Critically Ill Patients With COVID-19 [published online ahead of print, 2020 Mar 11]. *JAMA*. 2020
- <sup>6</sup> NHS. Clinical guide for the use of acute non-invasive ventilation in adult patients hospitalised with suspected or confirmed coronavirus during the coronavirus pandemic [19 March 2020, Version 1]. Retrieved from <https://www.england.nhs.uk/coronavirus/wp-content/uploads/sites/52/2020/03/clinical-guide-acute-niv-ventilation-v1-19-march-2020.pdf> March 22, 2020.
- <sup>7</sup> Handbook of COVID-19 Prevention and Treatment. The First Affiliated Hospital, Zhejiang University School of Medicine Compiled According to Clinical Experience. Retrieved from <https://covid-19.alibabacloud.com/> on March 18, 2020.
- <sup>8</sup> Xie, J., Tong, Z., Guan, X. et al. Critical care crisis and some recommendations during the COVID-19 epidemic in China. *Intensive Care Med* (2020). <https://doi.org/10.1007/s00134-020-05979-7>
- <sup>9</sup> Medicines & Healthcare Products Regulatory Agency (MHRA). Rapidly Manufactured Ventilator Systems (RMVS). Retrieved from [https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment\\_data/file/874279/RMVS001\\_Rapidly\\_Manufactured\\_Ventilator\\_Specification\\_PDF.pdf](https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/874279/RMVS001_Rapidly_Manufactured_Ventilator_Specification_PDF.pdf) March 23, 2020
- <sup>10</sup> US Department of Health and Human Services. Food and Drug Administration. Center for Devices and Radiological Health. Enforcement Policy for Ventilators and Accessories and Other Respiratory Devices During the Coronavirus Disease 2019 (COVID-19) Public Health Emergency. Guidance for Industry and Food and Drug Administration Staff. March 2020
- <sup>11</sup> Hick, JL, Hanfling D, Wynia MK, et al. Duty to Plan: Health Care, Crisis Standards of Care, and Novel Coronavirus SARS-CoV-2. *NAM Perspectives*. Discussion paper. National Academy of Medicine. Washington, DC. 2020. Retrieved from <https://doi.org/10.31478/202003b> March 17, 2020.
- <sup>12</sup> New York State Department of Health. Ventilator Allocation Guidelines. Retrieved from [https://www.health.ny.gov/regulations/task\\_force/reports\\_publications/docs/ventilator\\_guidelines.pdf](https://www.health.ny.gov/regulations/task_force/reports_publications/docs/ventilator_guidelines.pdf) March 17, 2020.
- <sup>13</sup> Phillips Healthcare. Respiratory Management for Severe COVID-19 Patients Webinar – Part1. Retrieved from <https://www.youtube.com/watch?reload=9&v=rMT97toZnJw&feature=youtu.be> March 23, 2020.
- <sup>14</sup> Simonds AK, Hanak A, Chatwin M, et al. Evaluation of droplet dispersion during non-invasive ventilation, oxygen therapy, nebuliser treatment and chest physiotherapy in clinical practice: implications for management of pandemic influenza and other airborne infections. *Health Technol Assess*. 2010
- <sup>15</sup> Cascella M, Rajnik M, Cuomo A, et al. Features, Evaluation and Treatment Coronavirus (COVID-19) [Updated 2020 Mar 8]. In: StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing; 2020 Jan-. Available from: <https://www.ncbi.nlm.nih.gov/books/NBK554776/>
- <sup>16</sup> Gonzalez J, Maisonneuve J, Oranger M, Mendoza-Ruiz A. Homecare respiratory equipment for patients suspected to be infected with the respiratory virus Covid19. Retrieved from <http://splf.fr/wp-content/uploads/2020/03/LES-PROCEDURES-DU-GAVO2-ProtectionVirale2020-MAJ14mars2020.pdf> March 18, 2020
- <sup>17</sup> CDC. Interim Infection Prevention and Control Recommendations for Patients with Suspected or Confirmed Coronavirus Disease 2019 (COVID-19) in Healthcare Settings. Retrieved from <https://www.cdc.gov/coronavirus/2019-ncov/infection-control/control-recommendations.html> March 17, 2020.
- <sup>18</sup> Lumis and AirCurve are essentially the same device, with different country variant naming, specific indications for use of variants of the Lumis and AirCurve device may vary by region.
- <sup>19</sup> Environmental Protection Agency. List N: disinfectants for use against SARS-CoV-2. EPA website. Retrieved from <https://www.epa.gov/pesticide-registration/list-n-disinfectants-use-against-sars-cov-2> March 22, 2020.
- <sup>20</sup> Australian Government Department of Health. Environmental cleaning and disinfecting principles – Version 1 (10/03/2020); Coronavirus disease (COVID-19). Retrieved from <https://www.health.gov.au/sites/default/files/documents/2020/03/environmental-cleaning-and-disinfection-principles-for-covid-19.pdf> Mar 18, 2020




**ResMed**

Rev 2.0

## APPENDIX – Device Specifications

Device availability will vary by region.


	<b>Astral 100/150</b>	<b>Stellar 100/150</b>	<b>GA ST-40P*</b>	<b>Lumis 150 ST-A</b>	<b>AirCurve ST-A</b>	<b>Lumis 150 ST</b>	<b>AirCurve ST</b>	<b>Flexo ST30*</b>
<b>Indications for use</b>	The Astral 100/150 devices provide continuous or intermittent ventilatory support for patients weighing more than 5 kg (11 lb) who require mechanical ventilation. The iVAPS mode with optional AutoEPAP is intended for patients weighing more than 30 kg (66 lb). The Astral device is intended to be used in home, institution/ hospital and portable applications for both invasive and non-invasive ventilation. The Astral device is not intended for use as an emergency transport ventilator.	The Stellar 100/150 is intended to provide ventilation for non-dependent, spontaneously breathing adult and paediatric patients (13kg/30 lb and above) with respiratory insufficiency, or respiratory failure, with or without obstructive sleep apnoea. The device is for noninvasive use, or invasive use (with the use of the ResMed Leak Valve).	The GA non-invasive ventilator is an assist ventilator and is intended to augment patient breathing. It is intended solely for spontaneously breathing individuals who require mechanical ventilation: patients with respiratory failure, chronic respiratory insufficiency, or obstructive sleep apnea in a hospital or other institutional settings under the direction of a physician. The ventilator is intended to support patients weighing 20 kg (44 lb) or greater	The Lumis 150 ST-A device is indicated to provide non-invasive ventilation for patients weighing more than 13 kg/30 lb or more than 30 kg/66 lb in iVAPS mode with respiratory insufficiency or obstructive sleep apnoea (OSA). It is intended for home and hospital use.	The AirCurve 10 ST-A is indicated to provide non-invasive ventilation for patients weighing more than 30 lb (13 kg) with respiratory insufficiency or obstructive sleep apnea (OSA). The AirCurve 10 ST-A is intended for home and hospital use.	The Lumis 150 VPAP ST device is indicated to provide non-invasive ventilation for patients weighing more than 13 kg/30 lb or more than 30 kg/66 lb in iVAPS mode with respiratory insufficiency or obstructive sleep apnoea (OSA). It is intended for home and hospital use.	The AirCurve 10 ST device is indicated for the treatment of obstructive sleep apnea (OSA) in patients weighing more than 66 lb (30 kg). It is intended for home and hospital use.	Flexo Bi-level positive airway pressure therapy device is used to provide non-invasive ventilatory support of mechanical ventilation for <ul style="list-style-type: none"> <li>• Patients who are independent and spontaneously breathing</li> <li>• Patients with respiratory insufficiency and sleep breathing disorders</li> </ul> It can provide both a stable continuous positive airway pressure and a bi-level positive airway pressure. It is not a life support ventilator. It is intended to be used in the home or professional medical environment.
<b>Device Type</b>	Mechanical ventilator for life support	Mechanical ventilator non-life support	Mechanical ventilator non-life support	Bi-level/BiPAP with alarms	Bi-level/BiPAP with alarms	Bi-level/BiPAP	Bi-level/BiPAP	Bi-level/BiPAP



**ResMed**

Rev 2.0


	<b>Astral 100/150</b>	<b>Stellar 100/150</b>	<b>GA ST-40P*</b>	<b>Lumis 150 ST-A</b>	<b>AirCurve ST-A</b>	<b>Lumis 150 ST</b>	<b>AirCurve ST</b>	<b>Flexo ST30*</b>
<b>Patient Weight</b>	>5 kg (all modes except iVAPS)	≥13 kg (all modes except iVAPS)	≥20 kg	>13 kg (all modes except iVAPS)	>13 kg (all modes except iVAPS)	>13kg (all modes except iVAPS)	>30kg	-
	>30 kg (iVAPS)	≥30 kg (iVAPS)		≥30 kg (iVAPS)	≥30 kg (iVAPS)	≥30kg iVAPS		
<b>Ventilation Delivery Specifications</b>								
<b>Patient Interface</b>	non-invasive invasive	non-invasive Invasive (with the use of ResMed leak valve)	non-invasive	non-invasive	non-invasive	non-invasive	non-invasive	non-invasive
<b>Circuit Type</b>	Single valve Single leak Double	single leak	single leak	single leak	single leak	single leak	single leak	single leak
<b>Pressure or Volume cycled</b>	Pressure Volume	Pressure	Pressure	Pressure	Pressure	Pressure	Pressure	Pressure
<b>Operating Pressure Range</b>	2-50 cmH <sub>2</sub> O: (S)T, P(A)C, iVAPS 3-20 cmH <sub>2</sub> O: CPAP 0-50 cmH <sub>2</sub> O: PS, P-SIMV, V-SIMV	2-40 cmH <sub>2</sub> O: S, ST, T, PAC, iVAPS 4-20 cmH <sub>2</sub> O: CPAP	4-40 cmH <sub>2</sub> O	2 - 30 cmH <sub>2</sub> O: S, ST, T, PAC, iVAPS 4-20 cmH <sub>2</sub> O: CPAP	3-30 cmH <sub>2</sub> O: S, ST, T, PAC, iVAPS 4-20 cmH <sub>2</sub> O: CPAP	2-25 cmH <sub>2</sub> O: S, ST, T, PAC, iVAPS 4-20 cmH <sub>2</sub> O: CPAP	3-25 cmH <sub>2</sub> O: S, ST, T 4-20 cmH <sub>2</sub> O: CPAP	4-30 cmH <sub>2</sub> O: S, T, ST, APCV 4-20 cmH <sub>2</sub> O: CPAP
<b>Therapy Modes</b>	CPAP, (S)T, P(A)C, (A)CV, P(A)CV, P-SIMV, V-SIMV, PS, iVAPS	CPAP, S, ST (optional iBR), T, PAC, iVAPS	CPAP, S, T, ST, APCV, CPAV, TVV-t, TVV-ST, TVV-APCV	CPAP, S, ST (optional iBR), T, PAC, iVAPS	CPAP, S, ST (optional iBR), T, PAC, iVAPS	CPAP, S, ST (optional iBR), T, PAC, iVAPS	CPAP, S, ST, T	CPAP, S, T, ST, APCV, TVV
<b>Performance Specifications</b>								
<b>Pressure Support Range</b>	-	0-38 cmH <sub>2</sub> O 0-30 cmH <sub>2</sub> O (in iVAPS mode)	0-38 cmH <sub>2</sub> O	0-28 cmH <sub>2</sub> O	0-27 cmH <sub>2</sub> O	0-23 cmH <sub>2</sub> O	0-22 cmH <sub>2</sub> O	0-26 cmH <sub>2</sub> O
<b>Intended Volume Range</b>	Adult: 100 - 2,500 mL, Paed: 50 - 300 mL: (A)CV, V-SIMV	50-3000 mL	200-2000 mL	100-2500 mL	100-2500 mL	100-2500 mL	100-2500 mL	50-2500 mL
<b>Supplementary Oxygen – Max Flow</b>	30 L/min	30 L/min	30 L/min	15 L/min (S, ST, T, PAC, CPAP) 4 L/min (iVAPS)	15 L/min (S, ST, T, PAC, CPAP) 4 L/min (iVAPS)	15 L/min (S, ST, T, PAC, CPAP) 4 L/min (iVAPS)	15 L/min (S, ST, T)	n/a



**ResMed**

Rev 2.0


	<b>Astral 100/150</b>	<b>Stellar 100/150</b>	<b>GA ST-40P*</b>	<b>Lumis 150 ST-A</b>	<b>AirCurve ST-A</b>	<b>Lumis 150 ST</b>	<b>AirCurve ST</b>	<b>Flexo ST30*</b>
<b>Supplementary Oxygen – where to add</b>	Device oxygen inlet	Device oxygen inlet	Device oxygen inlet	Device air outlet (eg. Side port of ClimateLineAir Oxy tube), mask	Device air outlet (eg. Side port of ClimateLineAir Oxy tube), mask	Device air outlet (eg. Side port of ClimateLineAir Oxy tube), mask	Device air outlet (eg. Side port of ClimateLineAir Oxy tube), mask	n/a
<b>SpO2 measurement</b>	Optional accessory	Optional accessory	Optional accessory	Optional accessory	Optional accessory	Optional accessory	Optional accessory	n/a
<b>FiO2 monitoring</b>	Optional accessory	Optional accessory	Optional accessory	n/a	n/a	n/a	n/a	n/a
<b>Internal Battery Runtime</b>	8 hours (a new battery under normal conditions)	3 hours (a new battery under normal conditions)	4 hours	n/a	n/a	n/a	n/a	n/a
	Lithium-Ion battery, 14.4 V, 6.6 Ah, 95 Wh	Lithium-Ion battery, 14.4 V, 2.75 Ah, 40 Wh						
<b>Humidification</b>	External	Integrated or External	External	Integrated or External	Integrated or External	Integrated or External	Integrated or External	External
<b>Safety</b>								
<b>Max single fault pressure</b>	90 cmH <sub>2</sub> O	60 cm H <sub>2</sub> O (in all modes)	-	30 cmH <sub>2</sub> O for >6 s, or	30 cmH <sub>2</sub> O for >6 s, or	30 cmH <sub>2</sub> O for >6 s, or	30 cmH <sub>2</sub> O for >6 s, or	-
				40 cmH <sub>2</sub> O for >1 s				
<b>Design Life</b>	Device: 8 years	Device: 5 years	Device: 10 years	Device, PSU: 5 years	Device, PSU: 5 years	Device, PSU: 5 years	Device, PSU: 5 years	Device: 5 years
		Air Tubing: 6 months		Cleanable Humidifier: 2.5 years				
				Air Tubing: 6 months				
<b>Alarms</b>	Total power failure	Power fail	High/low pressure	Power fail	Power fail	n/a	n/a	Patient circuit disconnect
	Circuit disconnection	Blocked tube	High/low oxygen pressure	Blocked tube	Blocked tube			
	Low Pressure, High Pressure	Tube disconnected, System fault	No oxygen supply	Tube disconnected	Tube disconnected			High inspiratory pressure
	Obstruction	High temperature, Internal battery low/empty	High/low inspiratory pressure	High leak	High leak			



**ResMed**

Rev 2.0

	<b>Astral 100/150</b>	<b>Stellar 100/150</b>	<b>GA ST-40P*</b>	<b>Lumis 150 ST-A</b>	<b>AirCurve ST-A</b>	<b>Lumis 150 ST</b>	<b>AirCurve ST</b>	<b>Flexo ST30*</b>
	Apnea	Over pressure	Patient circuit disconnection	Non-vented mask	Non-vented mask			Apnea
	Low / High MVe , Low / High MVi	High leak	Low/failed internal battery	Apnea	Apnea			
	Low / High Vte, Low / High Vti	Non-vented mask	Internal battery powering	Low SpO2 (when oximeter connected)	Low SpO2 (when oximeter connected)			Low tidal volume
	Low / High Resp rate	Low minute ventilation	external battery powering	Low minute ventilation	Low minute ventilation			
	High leak	Apnea	Power return	System fault	System fault			Low minute volume
	Ventilation stopped	High/Low pressure	Tube disconnected					
	Low / High SpO <sub>2</sub>	High/Low respiratory rate,	High extremity pressure					High/low leak
	Low / High pulse rate	High/Low FiO <sub>2</sub>	High/low external battery voltage					
	Low / High FiO <sub>2</sub>	Low SpO <sub>2</sub>	Pressure sensor failure					
	NV mask/Rebreathing		Communication error with main board					
	Incorrect circuit adapter		Cooling fan speed too low					
	Critically low battery		Blocked airway path					
	Incorrect circuit attached		High/low leak					
	Safety reset complete		Apnea					
	Battery inoperable		Low MV					
	Low / High PEEP		High/low BPM					
	Device overheating		High/low VT					
	Pressure line disconnected		High/low TVV					
	Self-test failed							
	Flow sensor not calibrated							



**ResMed**

Rev 2.0

	<b>Astral 100/150</b>	<b>Stellar 100/150</b>	<b>GA ST-40P*</b>	<b>Lumis 150 ST-A</b>	<b>AirCurve ST-A</b>	<b>Lumis 150 ST</b>	<b>AirCurve ST</b>	<b>Flexo ST30*</b>
No SpO2 monitoring								
No FiO2 monitoring								
Internal battery degraded								
Low internal battery								
Circuit fault								
Using internal battery								
Battery 1 fault, Battery 2 fault								
Power fault/No charging								
PEEP blower failure								
System Fault alarms								
<b>Environmental</b>								
<b>Operating Temperature</b>	32°F to 104°F (0°C to 40°C)	+32°F (0°C) to + 95°F (35°C)	+5°C to +40°C	+41°F to +95°F (+5°C to +35°C)	+5°C to +35°C			
<b>Operating Humidity</b>	5 to 93% non-condensing	10%–95% non-condensing	10 to 93% non-condensing	10 to 95% relative humidity, non-condensing	10%–95% non-condensing			
<b>Storage and Transport Temperature</b>	-13°F to 158°F (-25°C to 70°C)	-4°F (-20°C) to +140°F (60°C)	-20°C to +55°C	-4°F to +140°F (-20°C to +60°C)	-20°C to +55°C			
<b>Storage and Transport Humidity</b>	5 to 93% non-condensing	10%–95% non-condensing	10 to 93% non-condensing	5 to 95% relative humidity, non-condensing	5 to 95% relative humidity, non-condensing	5 to 95% relative humidity, non-condensing	5 to 95% relative humidity, non-condensing	10%–95% non-condensing

\*GA ST and Flexo ST are devices offered by Curative, a subsidiary of ResMed.